

DCMTK - Conformance #1139

Add minimal support for new DICONDE Thermography Storage SOP Classes

2024-09-23 17:32 - Jörg Riesmeier

<b>Status:</b>	Closed	<b>Start date:</b>	2024-09-23
<b>Priority:</b>	Normal	<b>Due date:</b>	
<b>Assignee:</b>	Jörg Riesmeier	<b>% Done:</b>	100%
<b>Category:</b>	Library and Apps	<b>Estimated time:</b>	1:00 hour
<b>Target version:</b>	3.6.9	<b>Compiler:</b>	
<b>Module:</b>	dcmdata, dcmnet		
<b>Operating System:</b>			
<b>Description</b>			
DICOM 2024d introduces two new Storage SOP Class UIDs for DICONDE-ASTM E3440-“Thermography”. Support for these UIDs and associated Storage SOP Classes should be added to the DCMTK.			

History

#1 - 2024-11-06 21:50 - Jörg Riesmeier

- Status changed from New to Closed
- Assignee set to Jörg Riesmeier
- % Done changed from 0 to 100
- Estimated time set to 1:00 h

Closed with commit 1a9d3c8d2.